

## ***Swing Aside Mount DSAM4***

The Swing Aside Mount combines the features of the original Locking Swing Away Mount with ROP components for use on a tilt seat. The mount can be swung away from the user by releasing the lock-pin connected to the cable that runs below the horizontal tube.

- For greater independence the user may release the lock pin to swing the mount away.
- Suitable for use on tilt seat systems.
- Non-directional.
- Quickly removed from the wheelchair.
- Includes standard TUSB Quick Release.



Item code: DSAM4-  
Introduced in 2008