

## Rigid Mount (DRM1)

The Rigid Mount provides a simple, adaptable and adjustable mounting solution. The mount consists of two stainless steel tubes joined by a rigid connection at a right angle.



### KEY FEATURES:

- Simple, sturdy and cost-effective
- Suitable for most standard wheelchairs
- Quickly and easily removed from the wheelchair
- Non-directional. Attaches to either side of the wheelchair with the same components
- Can be lifted and swung to the side
- Available in ROP version with lock receiver

### CAUTION (DRM1 & DRM1ROP):

The Rigid Mount can support equipment that is heavy enough to unbalance an unoccupied wheelchair when swung out to the side. It is recommended to remove the device and/or mount before the wheelchair is vacated.

- The DRM1 Rigid Mount is only suitable for limited Tilt Seat Systems. The vertical tube must not be installed or moved beyond 15° during tilt
- The DRM1ROP Rigid Mount ROP cannot be swung to the side. The mount must be unlocked and removed for access to the user or during transfers