

Mini Swing Mount (M75-SM)

The Mini Swing Mount is designed to give independence to users who can move a mount away unaided. The top mount tube rotates around 270 degrees with 7 stop points at each 45 degrees. Free of locking mechanisms or levers stop positions are friction held with a force of 2.9lbs/1.3kg at 14"/35cm required to push or pull the mount tube.



KEY FEATURES:

- The mount is non-directional and may be installed on either side of a wheelchair
- M75-SM standard side tube dimension 24" length, 6" offset, and horizontal tube 14" straight

DEFAULT CONFIGURATION:

- **Frame Clamp Inner Piece:**
1" standard UFC1000IP
- **Quick Release Base:**
MQR75
- **Frame Clamp Assembly Style:**
Side tube 24" length, 6" offset – M75SC-24x6

CAUTION:

- In the swung away position the mount and device may protrude beyond the side of the wheelchair making it vulnerable to collision if the chair is moved or driven.