


## Rolling Mount Offset Base Tall Height

DAESSY Rolling Mounts provide an independent floor stand solution for using an AAC/SGD communication device, tablet or similar device. Fitted with a DAESSY quick release base a Rolling Mount is an excellent companion to the wheelchair mounts offering easy transfer of the device from wheelchair to Rolling Mount for use in other locations such as a stander, bed, or recliner.

All Rolling Mounts have easy tool-free height adjustability and low-profile bases with double locking casters.

### ROM32T



- Inside clearance 30"
- Overall width 33 ¼" x overall depth 25 ½"
- Tall height range provides device height adjustment between 39" and 57" from the floor.

### ROM36T



- Inside clearance 34"
- Overall width 37 ¼" x overall depth 25 ½"
- Tall height range provides device height adjustment between 39" and 57" from the floor.

### KEY FEATURES:

- An offset base places the crossbar to the front to allow the base to fit around a wheelchair, positioner chair, or recliner, etc. Different base widths are available.
- The DAESSY TUSB quick release base is included by default, alternates available.

### CAUTION:

- The Rolling Mount Series must be used on a level floor with all wheels securely in contact with the floor before extending to full height.
- A limiter stops the horizontal tube of the mount to be swung outside of the Rolling Mount base for stability. Modifications or tampering with the limiter are not to be made by any person.